


EVEREST Titanium Ligating Clips

SAFE
SECURE
SMART



Everest Clip Simply Super ...

- Clip made of pure titanium.
- Transverse serrations guarantee no slipping from vessel.
- Chevron shape clip for complete vessel overlapping.
- MRI compatible & radio-opaque.
- Precise smart loading mechanism.
- Adhesive backing on the EVEREST cartridge.
- EVEREST Applicators with diamond dust jaws.

EVEREST Titanium Ligating Clips



Safety

EVEREST Clip made of pure high medical grade titanium ensure biocompatibility, malleability and compatibility with MRI.

Security

EVEREST Clip is secured on the ligated tissues by transverse and longitudinal serrations and the chevron shape of the clip guarantee complete tip to tip closure, security on the applicator's jaw ensured by our special diamond dust treatment of the applicator's grooves.

Smart

The smart cartridge design allows one hand frictionless loading with wide variety of color coded applicators to meet every need.

Adhesive cartridge backing holds on any surface.

EVEREST TITANIUM LIGATING CLIPS

Article	Clip size	Color	Packing
EVC010	Small	Red	30 cartridges of 6 clips
EVC020	Medium	Blue	30 cartridges of 6 clips
EVC030	Medium-Large	Green	20 cartridges of 6 clips
EVC040	Large	Orange	20 cartridges of 6 clips

EVEREST TITANIUM LIGATING CLIPS APPLIERS

Article	Clip size	Applier length	Description
EVC110	Small	15 cm - 6"	Curved Jaw
EVC111	Small	20 cm - 8"	Curved Jaw
EVC111A	Small	20 cm - 8"	Angled Jaw
EVC113	Small	28 cm - 11"	Curved Jaw
EVC010M	Small		Cardiac minimal invasive applier
EVC010L	Small		Laparoscopic 5 mm diameter applier
EVC120	Medium	15 cm - 6"	Curved Jaw
EVC121	Medium	20 cm - 8"	Curved Jaw
EVC121A	Medium	20 cm - 8"	Angled Jaw
EVC123	Medium	28 cm - 11"	Curved Jaw
EVC123A	Medium	28 cm - 11"	Angled Jaw
EVC020M	Medium		Cardiac minimal invasive applier
EVC020L	Medium		Laparoscopic 6.5 mm diameter applier
EVC121L	Medium		Laparoscopic 10 mm diameter applier
EVC131	Medium-Large	20 cm - 8"	Curved Jaw
EVC131A	Medium-Large	20 cm - 8"	Angled Jaw
EVC133	Medium-Large	28 cm - 11"	Curved Jaw
EVC133A	Medium-Large	28 cm - 11"	Angled Jaw
EVC030L	Medium-Large		Laparoscopic 10 mm diameter applier
EVC141	Large	20 cm - 8"	Curved Jaw
EVC141A	Large	20 cm - 8"	Angled Jaw
EVC143	Large	28 cm - 11"	Curved Jaw
EVC143A	Large	28 cm - 11"	Angled Jaw
EVC040L	Large		Laparoscopic 12.5 mm diameter applier

© 2014 MDD GmbH. All rights reserved. MDD GmbH reserves the right to change or modify specifications without notice.

MDD GmbH

Steigäcker 20 · D-78582 Balgheim · Germany

Tel.: +49(0)7424/92208 · Fax: +49(0)7424/92242

info@m-d-d.net · www.m-d-d.net

