



Product catalogue



2015


Contents

Drapes and sets	3
X-ray detectable towels, swabs and sponges	29
Minor procedure sets	33
Surgical gowns	45
Protective apparel	53
Surgical masks and caps	61
Gloves	69
Wound care	93
Bandages	101





Drapes and sets




Safety

Practicality


Comfort

- Safely strong – impermeable
- Plenty of absorbency – extra large absorbent patches
- Very low linting
- Very good drapeability
- Stays well in place, is not slippery
- Low reflections
- High puncture resistance – very durable

- 1 Large and very efficient absorbent patches, new structure
- 2 Blue/grey nonwoven, very low linting
- 3 A highly elastic multicomponent pe-film moisture barrier – puncture resistance
- 4 White backsheet with very low linting – more comfort to the patient, more strength



Complies to EN 13795.



Easily identifiable integrated tube holders in drapes with absorbent patch

Material chart



Absorbent and impermeable multilayer material



Absorbent and impermeable 2-layer material



Absorbent and impermeable reinforcement material (absorbent patch)



Polyethylene film



Elastic material




Adhesive edge



Universal set with Mayo stand cover, absorbent patches size Large

REF 1011-02


- 1 op-tape, 9 x 50 cm
- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 2 adhesive drapes 75 x 90 cm, absorbent patches Large 25 x 60 cm
- 1 adhesive drape medium 180 x 180 cm, absorbent patch Large 25 x 60 cm, tube holders
- 1 adhesive drape large 150 x 240 cm, absorbent patch Large 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/5/10



Universal set without Mayo stand cover, absorbent patches size Large

REF 1015-02


- 1 op-tape 9 x 50 cm
- 4 hand towels
- 2 adhesive drapes 75 x 90 cm, absorbent patches Large 25 x 60 cm
- 1 adhesive drape medium 180 x 180 cm, absorbent patch Large 25 x 60 cm, tube holders
- 1 adhesive drape large 150 x 240 cm, absorbent patch 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/6/12



Universal set with Mayo stand cover

REF 1005-02

- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 78 x 145 cm, reinforced, telescopically folded
- 2 hand towels
- 2 adhesive drapes 75 x 90 cm, 2-layer
- 1 adhesive drape medium 180 x 180 cm, 2-layer
- 1 adhesive drape large 150 x 240 cm, 2-layer
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/6/12




Drapes and sets

Universal sets

Universal set

REF 1006-01


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 op-tape 9 x 50 cm
- 2 adhesive drapes 75 x 90 cm, absorbent patch 15 x 50 cm
- 4 hand towels
- 1 adhesive drape Medium 180 x 180 cm, absorbent patch 15 x 50 cm, tube holders
- 1 adhesive drape Large 150 x 240 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/6/12



Paediatric universal set

REF 1030-02


- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 2 adhesive drapes 75 x 75 cm
- 1 adhesive drape medium 180 x 180 cm
- 1 adhesive drape large 150 x 240 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/7/14



Universal set with Mayo stand cover, absorbent patches size Medium

REF 1010-02


- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 2 adhesive drapes 75 x 90 cm, absorbent patches Medium 15 x 50 cm
- 1 adhesive drape medium 180 x 180 cm, absorbent patches Medium 15 x 50 cm, tube holders
- 1 adhesive drape large 150 x 240 cm
- 1 instrument table cover 150 x 190 cm
- Pack counts: 1/5/10



Split sheet set

REF 1050-02


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 1 op-tape 9 x 50 cm
- 1 adhesive drape 75 x 90 cm
- 1 adhesive split sheet 230 x 260 cm, split 20 x 100 cm, absorbent patch 75 x 140 cm, tube holders
- 1 adhesive drape large, 150 x 240 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/4/8



Split sheet set

REF 1050-50


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 op-tape 9 x 50 cm
- 4 hand towels
- 1 adhesive drape sheet 75 x 90 cm
- 1 adhesive split sheet 230 x 260 cm, split 20 x 85 cm
- 1 adhesive drape Large 150 x 240 cm
- 1 instrument table cover, size 150x190 cm
- Pack count: 1/5/10



Split sheet set

REF 1051-02

- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 1 op-tape 9 x 50 cm
- 1 adhesive drape 75 x 90 cm
- 1 adhesive split sheet 230 x 260 cm, split 20 x 100 cm, absorbent patch 75 x 140 cm, tube holders, split sheet folded so that the bottom of the split is visible
- 1 adhesive drape large, 150 x 240 cm, absorbent patch 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/4/8




Drapes and sets

Orthopaedic sets

Hip set, Standard

REF 1210-02


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 1 adhesive drape medium 180 x 180 cm, absorbent patch 25 x 60 cm, tube holders
- 1 elastic stockinette 30 x 120 cm
- 1 elastic support bandage 12 cm x 6 m
- 2 adhesive drapes 75 x 90 cm, absorbent patch 25 x 60 cm
- 3 op-tapes 9 x 50 cm
- 1 adhesive split sheet 230 x 260 cm, split 20 x 100 cm, absorbent patch 75 x 140 cm, tube holders
- 1 adhesive drape 150 x 240 cm, absorbent patch 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/3/6



HIP set, Large

REF 1215-02


- 2 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 1 adhesive drape medium 180 x 180 cm, absorbent patch 25 x 60 cm, tube holders
- 1 elastic stockinette 30 x 120 cm
- 1 elastic support bandage 12 cm x 6 m
- 1 instrument table cover 100 x 150 cm, folded
- 2 utility drapes 75 x 80 cm
- 3 op-tapes 9 x 50 cm
- 3 adhesive drapes 75 x 90 cm, absorbent patch 25 x 60 cm
- 1 fluid collection pouch 30 x 40 cm
- 1 adhesive split sheet 230 x 260 cm, split 20 x 100 cm, absorbent patch 75 x 140 cm, tube holders
- 1 adhesive drape large, 150 x 240 cm, absorbent patch 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/3/6



Hip set, Special

REF 1220-02


- 2 Mayo stand covers 78 x 145 cm, reinforced
- 4 hand towels
- 1 adhesive drape medium 180 x 180 cm, absorbent patch 25 x 60 cm, tube holders
- 1 elastic support bandage 12 cm x 6 m
- 1 elastic stockinette 30 x 120 cm
- 5 op-tapes 9 x 50 cm
- 3 adhesive drapes 75 x 90 cm, absorbent patch 25 x 60 cm
- 1 adhesive split sheet, 230 x 260 cm, split 20 x 100 cm, abs. patch 75 x 140 cm, tube holders
- 1 adhesive drape large, 150 x 240 cm, absorbent patch 25 x 60 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/3/6



Hip set

REF 1205-01


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 utility drape 100 x 150 cm
- 2 op-tape 9 x 50 cm
- 1 stockinette, 36 x 120 cm
- 1 adhesive split sheet 230 x 260 cm, absorbent patch 50 x 100 cm, split 20 x 85 cm
- 1 adhesive drape 150 x 300 cm, absorbent patch 15 x 50 cm, tube holders
- 1 Instrument table cover 150 x 190 cm
- Pack count: 1/4/8



Vertical drape set, Small

REF 1225-02


- 4 hand towels
- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 vertical drape 330 x 240 cm, incision film 40 x 30 cm, integrated fluid collection pouch with sieve and bottom outlet, 2 tube holders and 2 two-compartment suction/diathermy pouches, 24 x 24 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/7/14



Vertical drape set, Large

REF 1226-02


- 4 hand towels
- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 vertical drape 330 x 240 cm, incision film 86 x 25 cm, integrated fluid collection pouch with sieve and bottom outlet, 2 tube holders and 2 two-compartment suction/diathermy pouches, 24 x 24 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/7/14



Extremity set

REF 1230-02


- 4 hand towels
- 1 Mayo stand cover 78 x 145 cm, reinforced
- 2 op-tapes 9 x 50 cm
- 1 instrument table cover 150 x 190 cm, folded
- 1 adhesive drape sheet 75 x 90 cm
- 1 extremity drape 230 x 300 cm, elastic aperture Ø 7 cm, absorbent patch 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/4/8



Extremity set

REF 1231-02


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 1 instrument table cover 150 x 190 cm, folded
- 2 op-tapes 9 x 50 cm
- 1 elastic stockinette 22 x 75 cm
- 1 extremity drape 230 x 300 cm, elastic aperture Ø 7 cm, absorbent patch 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/4/8



Extremity set

REF 1231-100


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 op-tape 9 x 50 cm
- 4 hand towels
- 1 stockinette 36 x 65 cm
- 1 extremity drape 230 x 300 cm, absorbent patch 75 x 80 cm, tube holders, elastic aperture Ø 7 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/6/12



Extremity set

REF 1233-02

- 1 Mayo stand cover, 78 x 145 cm, reinforced
- 2 op-tapes 9 x 50 cm
- 4 hand towels
- 1 instrument table cover 150 x 190 cm, folded
- 1 adhesive drape 75 x 90 cm
- 1 elastic support bandage 10 cm x 4,5 m
- 1 extremity drape 230 x 300 cm, elastic aperture Ø 7 cm, absorbent patch 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/4/8




Drapes and sets

Orthopaedic sets

Hand/foot set

REF 1232-02


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 extremity drape 230 x 300 cm, elastic aperture \varnothing 3,5 cm, absorbent patch 50 x 100 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/6/12



Hand/foot set

REF 1232-50


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 op-tape 9 x 50 cm
- 4 hand towels
- 1 hand/foot drape 230 x 300 cm, elastic aperture \varnothing 3,5 cm, absorbent patch 75 x 80 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/6/12



Knee arthroscopy set

REF 1235-02


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 2 op-tapes 9 x 50 cm
- 1 instrument table cover 150 x 190 cm, folded
- 1 elastic stockinette 22 x 75 cm
- 1 knee arthroscopy drape with fluid collection pouch, 230 x 320 cm, elastic apertures \varnothing 5 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/4/8



Knee arthroscopy set

REF 1235-50


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 tube cover 15 x 240 cm, perforated tip
- 2 hand towels
- 1 op-tape 9 x 50 cm
- 1 stockinette, 36 x 65 cm
- 1 extremity drape 230 x 300 cm, elastic aperture \varnothing 7 cm, absorbent patch 50 x 100, tube holders
- 1 instrument table cover 150 x 140 cm
- Pack count: 1/5/10



Extremity set, bilateral

REF 1234-01

- 4 hand towels
- 2 op-tapes 9 x 50 cm
- 1 extremity drape 200 x 300 cm, bilateral, elastic apertures \varnothing 6 cm, absorbent patch 90 x 150 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/8/24




Urology and gynaecology sets

Gyn/Cyst set

REF 1250-01


- 1 utility drape, 75 x 80 cm
- 2 leggings, 75 x 120 cm
- 1 gyn/cyst drape 75 x 175 cm, aperture 9 x 15 cm
- 1 instrument table cover, 150 x 190 cm
- Pack count: 1/5/20



Gyn/Cyst set

REF 1250-50


- 2 leggings, 75 x 120 cm
- 1 gyn/cyst drape, 100 x 175 cm, adhesive aperture 9 x 15 cm
- 1 instrument table cover, 100 x 150 cm
- Pack count: 1/5/20



Gynaecology set

REF 1251-02


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 2 hand towels
- 1 tube holder 2,5 x 30 cm
- 1 gynaecology drape with integrated fluid collection pouch, 60 x 120 cm, aperture 9 x 12 cm.
- 1 gynaecology drape 270/230 x 260 cm, aperture 24 x 21 cm, integrated leggings
- 1 Instrument table cover 150 x 190 cm
- Pack count: 1/6/12



Gynaecology set

REF 1253-01

- 3 hand towel
- 1 drape sheet 75 x 120 cm
- 2 leggings 75 x 120 cm
- 1 gyn upper drape 115 x 150 cm
- 1 gyn lower drape 140 x 150 cm
- 1 instrument table cover 100 x 150 cm
- Pack count: 1/4/16




Urology and gynaecology sets

Abdominal-perineal set

REF 1248-50


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 2 op-tape 9 x 50 cm
- 2 hand towels
- 1 gyn-abdominal drape 220/240 x 335 cm, all in-one with integrated leggings and fluid collection pouch, with 2 adhesive apertures: top aperture 25 x 30, bottom aperture 9 x 15 cm (with integrated drape to cover the aperture when not needed)
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/5/10



Delivery set

REF 1255-01


- 1 instrument table cover (for under the mother), 90 x 150 cm
- 1 drape (stitching sheet) 75 x 120 cm
- 1 bed pad 60 x 90 cm
- 1 absorbent sheet (baby sheet) 56 x 75 cm
- 6 towels 21 x 25 cm
- 1 bed pad 60 x 60 cm
- Pack count: 1/12



C-section set

REF 1260-02

- 4 hand towels
- 1 op-tape 9 x 50 cm
- 1 Mayo stand cover 78 x 145 cm, reinforced
- 1 baby sheet 75 x 120 cm
- 1 C-section drape 186/250 x 330 cm, fluid collection pouch, incision film, Ø 16 x 18 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/5/10




Urology and gynaecology sets

C-section set

REF 1261-02


- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 1 op-tapes 9 x 50 cm
- 1 baby sheet 75 x 120 cm
- 1 C-section drape 330 x 230 cm, fluid collection pouch, integrated incision film 36 x 36 cm and tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/5/10



TUR set

REF 1270-02


- 1 op-tape 9 x 50 cm
- 1 utility drape 75 x 80 cm
- 2 leggings 75 x 120 cm
- 1 TUR-drape 75 x 175 cm
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/4/16



TUR set

REF 1272-01


- 1 tube holder velcron 2,5 x 30 cm
- 4 hand towels
- 1 op-tape 9 x 50 cm
- 1 TUR-drape 175/270 x 170 cm all-in-one, with integrated finger cot, aperture 6 and self-adhesive aperture 8 cm, integrated fluid collection pouch with sieve and outlet
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/6/12



ENT set

REF 1238-02


- 4 hand towels
- 1 Mayo stand cover 78 x 145 cm, reinforced
- 2 op-tapes 9 x 50 cm
- 1 adhesive drape 75 x 90 cm, absorbent patch 15 x 50 cm, tube holders
- 1 utility drape 75 x 80 cm
- 1 adhesive split sheet 150 x 250 cm split 7 x 70 cm, absorbent patch 50 x 50 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/5/10



Angiography set

REF 1240-02


- 1 instrument table cover 100 x 150 cm, folded
- 2 hand towels
- 1 absorbent sheet 50 x 80 cm
- 1 adhesive drape 50 x 60 cm
- 1 angiography drape 230 x 370 cm, with transparent side panels, absorbent patch 75 x 70 cm
- 1 instrument table cover 150 x 190 cm
- Pack count:1/5/10



Angiography set

REF 1241-01


- 1 instrument table cover 150 x 140 cm, folded
- 2 hand towels
- 1 absorbent drape 50 x 80 cm
- 1 adhesive drape 50 x 60 cm
- 1 angiography drape, 230 x 330 cm, with transparent side panels, apertures 5 x 7 cm (2 pcs) and 7 x 10 cm (2 pcs)
- 1 instrument table cover 150 x 140 cm
- Pack count:1/5/10



Dental set

REF 1288-01


- 1 dental drape 150 x 150 cm, diamond aperture 9 x 15 cm
- 2 op-tape 9 x 25 cm, blue
- 2 tube covers 8 x 120 cm
- 2 utility drapes 45 x 75 cm
- 1 instrument table cover 75 x 100 cm
- Pack count: 1/9/36



Dental set

REF 1290-02


- 1 dental drape 100 x 150 cm, aperture 9 x 13 cm
- 1 op-tape 9 x 50 cm
- 2 tube covers 8 x 120 cm
- 1 instrument table cover 75 x 100 cm
- Pack count: 1/14/56



Head/Neck drape

REF 1291-02


- 2 head/neck drapes 175 x 127 x 150 cm and 103 x 106 cm
- 1 instrument table cover 78 x 70 cm
- 1 op-tape 9 x 50 cm
- Pack count: 1/10/40



Laparoscopy set

REF 1267-02

- 1 Mayo stand cover 78 x 145 cm, reinforced
- 4 hand towels
- 2 op-tapes 9 x 50 cm
- 1 laparoscopy drape 230 x 300 cm, aperture 31 x 35 cm, absorbent patch 75 x 80 cm, tube holders
- 1 instrument table cover 150 x 190 cm
- Pack count: 1/7/14




Drapes and sets

Ophthalmic sets

Ophthalmic set

REF 1280-01


- 1 hand towel
- 1 ophthalmic drape 150 x 140 cm, with fluid collection pouch and integrated incision film
- 1 instrument table cover 150 x 140 cm
- Pack count: 1/7/28



Ophthalmic set

REF 1281-01


- 1 hand towel
- 1 ophthalmic drape 150 x 140 cm, fluid collection pouch and integrated incision film with oval aperture 3,5 x 6,5 cm
- 1 instrument table cover 150 x 140 cm
- Pack count: 1/7/28



Ophthalmic set

REF 1282-01


- 1 incision film 10 x 12 cm
- 1 ophthalmic drape 150 x 180 cm, fluid collection pouches, oval aperture 5 x 7 cm
- 1 instrument table cover 150 x 140 cm
- Pack count: 1/7/28




Adhesive split sheet

REF 1308

- 150 x 183 cm, split 7 x 70 cm
- Pack count: 1/9/36



150 x 183 cm
U = 7 x 70 cm



230 x 290 cm
U = 6 x 90 cm

Adhesive split sheet bilateral


REF 1309

- 230 x 290 cm, splits 6 x 90 cm absorbent patch 150 x 100 cm, tube holders
- Pack count: 1/6/12


Adhesive split sheets

REF 1310-01

- 150 x 250 cm, split 7 x 70 cm, absorbent patch 50 x 50 cm, tube holders
- Pack count: 1/10/40



150 x 250 cm
U = 7 x 70 cm




230 x 260 cm
U = 20 x 100 cm

REF 1315-01


- 230 x 260 cm, split 20 x 100 cm, absorbent patch 75 x 140 cm, tube holders
- Pack count: 1/8/16

REF 1316-01

- 230 x 260 cm, split 20 x 100 cm, absorbent patch 75 x 140 cm, tube holders, folded so that bottom of the split is visible
- Pack count: 1/8/16



230 x 260 cm
U = 20 x 100 cm



200 x 260 cm
U = 7 x 100 cm


REF 1318-01

- 200 x 260 cm, split 7 x 100 cm, absorbent patch 75 x 140 cm
- Pack count: 1/9/18

Hand/foot drapes

REF 1320-01


- 150 x 180 cm, elastic aperture Ø 3,5 cm, absorbent patch 50 x 100 cm, tube holders
- Pack count: 1/17/34



150 x 180 cm

REF 1324-01

- 230 x 300 cm, elastic aperture Ø 3,5 cm, absorbent patch 50 x 100 cm, tube holders
- Pack count: 1/8/16



230 x 300 cm


Drapes and sets

Speciality drapes

Knee arthroscopy drape

REF 1328


- 230 x 320 cm, with fluid collection pouch, elastic apertures \varnothing 5 cm
- Pack count: 1/7/14



Extremity drape

REF 1325-01


- 230 x 300 cm, elastic aperture \varnothing 7 cm, absorbent patch 50 x 100 cm, tube holders
- Pack count: 1/8/16



Vertical drapes


REF 1330-02

- 330 x 240 cm, incision film 40 x 30 cm, integrated fluid collection pouch with sieve and bottom outlet, 2 tube holders and 2 two-compartment suction/diathermy pouches, 24 x 24 cm
- Pack count: 1/14/28



REF 1332-02


- 330 x 240 cm, incision film 86 x 25 cm, integrated fluid collection pouch with sieve and bottom outlet, 2 tube holders and 2 two-compartment suction/diathermy pouches, 24 x 24 cm
- Pack count: 1/10/20



Lower vascular drape (Split sheet with genital cover)

REF 1317-01

- 230 x 260 cm split 20 x 100 cm
- Pack count: 1/12/24




230 x 260 cm
U = 20 x 100 cm

Angiography drape

REF 1340-01

- 230 x 300 cm, transparent side panels, absorbent patch 75 x 70 cm
- Pack count: 1/15/30




230 x 300 cm

Perineal drape

REF 1356-01

- 75 x 120 cm, aperture 8 x 12 cm, splits and adhesive tapes at edges
- 75 x 90 cm utility drape
- Pack count: 1/17/68



Laparoscopy drapes

REF 1365-01


- 230 x 300 cm, triangular aperture 25 x 20 cm
- Pack count: 1/9/18

REF 1366-01

- 230 x 300 cm, aperture 20 x 31 cm
- Pack count: 1/9/18



230 x 300 cm




230 x 300 cm

Thoracic drape

REF 1375-01

- 220/300 x 330 cm, aperture 32 x 40 cm, with incision film and fluid collection pouch
- Pack count: 1/7/14



200/300 x 300 cm

Under buttock drape

REF 1741-01

- 75 x 100 cm
- Pack counts: 1/25/100




75 x 100 cm

Dislocation bag

REF 1878

- 120 x 75 cm
- Pack counts: 1/20/80




120 x 75 cm

Utility drapes and adhesive drapes

Utility drapes

REF 1408


45 x 75 cm
Pack count: 1/30/120



45 x 75 cm



50 x 60 cm



75 x 90 cm

REF 1410-01

50 x 60 cm
Pack count: 1/45/180


REF 1415-01

75 x 90 cm
Pack count: 1/25/100


Adhesive drape sheets

REF 1456-01

37 x 45 cm
Pack count: 1/60/240



37 cm




50 cm

REF 1452-01


50 x 50 cm
Pack count: 1/30/120

REF 1455-01

50 x 60 cm
Pack count: 1/40/160



50 x 60 cm




45 cm

REF 1451-01


45 x 75 cm
Pack count: 1/30/120

REF 1461-01

75 x 75 cm
Pack count: 1/15/60



75 cm




75 x 90 cm

REF 1460-02


75 x 90 cm
Pack count: 1/16/64

REF 1454-01

75 x 100 cm
Pack count: 1/15/60



75 cm



100 cm

REF 1453-01


100 x 100 cm
Pack count: 1/13/52

Adhesive, aperture and adhesive aperture drapes

Absorbent drape

REF 1463-02


70 x 60 cm
Pack count: 1/18/72



Paediatric absorbent drape

REF 1463-50


70 x 60 cm
Pack count: 1/40/160



Adhesive drape medium

REF 1465-02


180 x 180 cm
Pack count: 1/8/32




Adhesive drapes large

REF 1464

180 x 254 cm
Pack count: 1/7/28



180 x 254 cm



150 x 240 cm


REF 1470-02

150 x 240 cm
Pack count: 1/8/32

Aperture drapes

REF 1555-01


50 x 60 cm, Ø 6 x 8 cm
Pack count: 1/45/180



50 x 60 cm
Ø 6 x 8 cm

REF 1565-01

75 x 90 cm, Ø 6 x 8 cm
Pack count: 1/25/100




75 x 90 cm
Ø 6 x 8 cm

Adhesive aperture drapes

REF 1505-01


50 x 60 cm, Ø 6 x 8 cm
Pack count: 1/45/180



50 x 60 cm
Ø 6 x 8 cm

REF 1510-01

75 x 80 cm, Ø 6 x 15 cm (epidural drape)
Pack count: 1/25/100




75 x 80 cm
Ø 6 x 15 cm

Adhesive aperture drape transparent

REF 1516

80 x 130 cm, Ø 15 x 6 cm
Pack count: 1/35/140




80 x 130 cm
Ø 15 x 6 cm

Aperture drapes and adhesive aperture drapes


Adhesive aperture drapes

REF 1515-01


75 x 90 cm, Ø 6 x 8 cm
Pack count: 1/25/100



75 x 90 cm
Ø 6 x 8 cm



100 x 120 cm
Ø 11 cm



120 x 150 cm
Ø 5 x 7 cm

REF 1517


100 x 125 cm, Ø 11cm
Pack count: 1/20/80

REF 1520-01

120 x 150 cm, Ø 5 x 7 cm
Pack count: 1/20/80


REF 1525-01

150 x 180 cm, Ø 5 x 7 cm
Pack count: 1/14/56



REF 1530-01

150 x 180 cm, Ø 6 x 15 cm
Pack count: 1/12/48




REF 1535-02

230 x 250 cm, Ø 12 cm
Pack count: 1/7/28

150 x 180 cm
Ø 5 x 7 cm


150 x 180 cm
Ø 6 x 15 cm



230 x 250 cm
Ø 12 cm


REF 1537-01

200 x 280 cm, Ø 14 x 14 cm
Pack count: 1/7/28



REF 1538-01

200 x 240 cm, Ø 12 cm
Pack count: 1/9/36



REF 1539-01

200 x 280 cm, Ø 15 cm
Pack count: 1/7/28

200 x 280 cm
Ø 14 x 14 cm


200 x 240 cm
Ø 12 cm

200 x 280 cm
Ø 15 cm

Variable aperture drapes


1512-01

50 x 75 cm
Pack count: 1/35/140



1513-01

75 x 100 cm
Pack count: 1/20/80



Paediatric adhesive aperture drape, variable

1511

37 x 65 cm
Pack count: 1/25/100





Table covers and components

Instrument table covers

REF 1715-01


100 x 150 cm
Pack count: 1/18/72



100 x 150 cm

REF 1720-01


140 x 150 cm
Pack count: 1/16/64




140 x 150 cm

REF 1725-01

150 x 190 cm
Pack count: 1/15/60



150 x 190 cm




150 x 190 cm


Instrument table cover, reinforced

REF 1727-01

150 x 190 cm
Pack count: 1/10/40



78 x 145 cm



78 x 145 cm

Mayo stand cover, reinforced

REF 1712-01

78 x 145 cm
Pack count: 1/20/80

Mayo stand cover, extra reinforced

REF 1714-01

78 x 145 cm
Pack count: 1/20/80

Adhesive Op-tapes

REF 1495-02

9 x 50 cm
Pack count: 1/100/400



9 x 50 cm

REF 1496-02

9 x 50 cm, 2 pcs/bag
Pack count: 2/120/480



9 x 50 cm

2 pcs

Tube holder

REF 1498-01

2,5 x 30 cm



2,5 x 30 cm

Seat cover

REF 1730-01

Ø 90 cm
Pack count: 1/7/28





Table covers and components

Pack count: 1/60/240

Arm cover

REF 1398-01

Pack count: 1/30/120




Leggings

REF 1390-01

75 x 120 cm, 2 pcs/bag

Pack count: 2/18/72


75 x 120 cm
2 pcs

Genital cover

REF 1393-01

20 x 42 cm

Pack count: 1/55/220



20 x 42 cm

Elastic stockinettes

REF 1395-01

30 x 120 cm, size L

Pack count: 1/12/48



30 x 120 cm



22 x 75 cm

REF 1396

22 x 75 cm, size S


Pack count: 1/21/84

Fluid collection pouch

REF 1870

30 x 40 cm

Pack count: 1/55/220




30 x 40 cm

Suction and diathermy bag

REF 1872

30 x 40 cm

Pack count: 1/55/220




30 x 40 cm

Fluid collection pouch

REF 1873

62 x 50 cm, sieve and outlet

Pack count: 1/11/44





62 x 50 cm

Equipment drapes

Tube covers

REF 1810

8 x 120 cm, simple tip
Pack count: 1/65/260





REF 1815

15 x 250 cm, simple tip
Pack count: 1/30/120


REF 1820

15 x 240 cm,
Perforated tip
Pack count: 1/30/120



REF 1825

17 x 240 cm,
Elastic tip
Pack count: 1/30/120




REF 1830

15 x 230 cm, ring insert
Narrow necked simple tip
Pack count: 1/12/48

Banded bags


REF 1850

Ø 90 cm
Pack count: 1/18/72



REF 1855-01

Ø 140 cm
Pack count: 1/10/40



Label information in accordance with EN ISO 15223-1 and EN 1041

evercare® is a registered trademark of OneMed Group Oy

Article number

In all products four (4) removable labels for traceability


Sterilization method:
– ethylene oxide EO

Single use

CE label: certifies that the product complies with the European legal requirements for medical devices

The identification number of the notified body

Manufacturers contact information



LOT code in all (4) removable labels

Use before

Barcodes in two (2) removable labels

Size

Does not contain natural rubber latex