



EVEREST

Titanium Ligating Clips





AARUT™
Quality.Conviction.Vision

EVEREST Titanium Ligating Clips

safe | secure | smart | simply super



EVEREST CLIP simply super Advantages at a glance

- Clip made of pure titanium
- Transverse serrations guarantee no slipping from vessel
- Chevron shape clip for complete vessel overlapping
- MRI compatible & radio-opaque
- Precise smart loading mechanism
- Adhesive backing on the EVEREST cartridge
- EVEREST appliers with diamond dust jaws



EVEREST Titanium Ligating Clips

safe | secure | smart



safe

EVEREST Clip made of pure high medical grade titanium ensure biocompatibility, malleability and compatibility with MRI.

secure

EVEREST Clip is secured on the ligated tissues by transverse and longitudinal serrations and the chevron shape of the clip guarantee complete tip to tip closure, security on the applier´s jaw ensured by our special diamond dust treatment of the applier´s grooves.

smart

The smart cartridge design allows one hand frictionless loading with wide variety of color coded appliers to meet every need. Adhesive cartridge backing holds on any surface.

Article	Clip size	Color	Packing
EVC060	micro	white	30 cartridges of 6 clips
EVC010	small	red	30 cartridges of 6 clips
EVC010C	small	red	multipack 15 x 24 clips
EVC011	small	yellow	30 cartridges of 6 clips
EVC010C-Y	small	yellow	multipack 15 x 24 clips
EVC050	small-medium	purple	30 cartridges of 6 clips
EVC020	medium	blue	30 cartridges of 6 clips
EVC020C	medium	blue	multipack 15 x 24 clips
EVC030	medium-large	green	20 cartridges of 6 clips
EVC040	large	orange	20 cartridges of 6 clips

Article	Clip size	Applier length	Description
EVC160	micro	15 cm	curved jaws
EVC161	micro	20 cm	curved jaws
EVC163	micro	28 cm	curved jaws
EVC110	small	15 cm	curved jaws
EVC111	small	20 cm	curved jaws
EVC111A	small	20 cm	angled jaws
EVC113	small	28 cm	curved jaws
EVC113A	small	28 cm	angled jaws
EVC010L	small	Laparoscopic applier 5 mm diameter	
EVC010M	small	Minimal invasive applier	
EVC110Y	small	15 cm	curved jaws
EVC111-Y	small	20 cm	curved jaws
EVC111A-Y	small	20 cm	angled jaws
EVC113-Y	small	28 cm	curved jaws
EVC113A-Y	small	28 cm	angled jaws
EVC150	small-medium	15 cm	curved jaws
EVC151	small-medium	20 cm	curved jaws
EVC120	medium	15 cm	curved jaws
EVC121	medium	20 cm	curved jaws
EVC121A	medium	20 cm	angled jaws
EVC123	medium	28 cm	curved jaws
EVC123A	medium	28 cm	angled jaws
EVC020L	medium	Laparoscopic applier 6.5 mm diameter	
EVC121L	medium	Laparoscopic applier 10 mm diameter	
EVC020M	medium	Minimal invasive applier	
EVC131	medium-large	20 cm	curved jaws
EVC131A	medium-large	20 cm	angled jaws
EVC133	medium-large	28 cm	curved jaws
EVC133A	medium-large	28 cm	angled jaws
EVC030L	medium-large	Laparoscopic applier 10 mm diameter	
EVC141	large	20 cm	curved jaws
EVC141A	large	20 cm	angled jaws
EVC143	large	28 cm	curved jaws
EVC143A	large	28 cm	angled jaws
EVC040L	large	Laparoscopic applier 12.5 mm diameter	

ABARUT GmbH & Co. KG



ABARUT™
Quality.Conviction.Vision

ABARUT GmbH & Co. KG
Pforzheimer Str.128 76275 Ettlingen GERMANY
Web: www.abarut.de Email : info@abarut.de